



Delta with Electric Adjustment

Enabling carer or user, where appropriate, to adjust the tilt-in-space facility and/or the elevating legrest by the simple touch of a button provides an electrically powered motion to both these mechanisms.

Utilising an undemanding hand controller the movement of the chair's tilt-in-space and legrest is smooth and gentle. The Delta with Electric adjustment reduces effort whilst increasing the positioning range and flexibility of the use of the chair.

Delta with Gas Action Adjustment

Effective and effortless movement for the carer, comfort and support for the user is key to the Delta range.

The Delta with Gas action adjustment enables fluid and precise movement of the tilt-in-space facility by the release of a lever situated on the push handle.

The tilt-in-space gently moves back and forth as directed by the carer with minimal effort whilst the legrest can be elevated by a separate release mechanism at the front of the chair. The legrest gently raises or lowers accordingly in a smooth arcing motion without affecting the positioning of the footboard.